

Protect From Light Chart | Laboratory Specimen

Mayo Test ID	Test Name
VITC	Ascorbic Acid, Plasma
BCARO	Beta-Carotene, Serum
AFBIL	Bilirubin, Amniotic Fluid
CDP	Chlordiazepoxide, Serum
CI	Inherited Conjugated Hyperbilirubinemias, Urine
MTHX	Methotrexate, Serum
PBGU	Porphobilinogen, Quantitative, Random, Urine
PTP	Porphyryns, Profile, Plasma
PQNRU	Porphyryns, Quantitative, Random, Urine
PQNU	Porphyryns, Quantitative, Urine
PLP	Pyridoxal 5-Phosphate (PLP), Plasma
VITB2	Riboflavin (Vitamin B2), Plasma
B6PRO	Vitamin B6 Profile (PLP and PA), Plasma
VITE	Vitamin E, Serum (9/5/2018)

Each specimen needs to be protected from light at all times. Please securely wrap the entire specimen in aluminum foil or paper towel or any other light protecting material.